

Laboratory 6 Sample Protocol

Contents

- Momentum Transfer in a Jet
- Objectives
- Safety considerations
- Procedure
- Data Acquisition
- Data Analysis
- Discussion/Interpretation
- Data Records

Momentum Transfer in a Jet

Objectives

- To demonstrate that the force on a vane is proportional to the rate of delivery of momentum
- To show that you can predict the force on a vane from a combination of its surface shape and the properties of the jet directed at it

Safety considerations

1. Don't drink the water!
2. Don't come to lab intoxicated!
3. Don't wear shorts, skirt.

4. Don't wear sandals; shoes or boots (closed toe, Crocks don't count)
5. What about electricity and water?
6. What about moving water and your eyeballs?

Procedure

Flat Plate Vane

0. Install the Vane
1. Ensure the weight beam is balanced at the zero position
2. Start the hydraulic bench and adjust flow (work from high flow to low flow)
3. Move the jockey weight until the beam balances again
 - Record distance, y , from the zero position
 - Record the flow rate using the hydraulic bench
4. Change the hydraulic bench flow rate in relatively equal increments for 4 additional readings
 - Record distance, y , from the zero position for each trial
 - Record the flow rate using the hydraulic bench for each trial

Hemispherical Cup Vane

0. Install the Cup
1. Ensure the weight beam is balanced at the zero position
2. Start the hydraulic bench and adjust flow (work from high flow to low flow)
3. Move the jockey weight until the beam balances again
 - Record distance, y , from the zero position
 - Record the flow rate using the hydraulic bench
4. Change the hydraulic bench flow rate in relatively equal increments for 4 additional readings
 - Record distance, y , from the zero position for each trial
 - Record the flow rate using the hydraulic bench for each trial

Data Acquisition

Follow the procedure above, be sure to record your trails in the tables below

Flat Plate

Trial	Distance (y)	Volume (ml)	Time (s)
1			
2			
3			
4			
5			

Hemispherical Cup

Trial	Distance (y)	Volume (ml)	Time (s)
1			
2			
3			
4			
5			

Data Analysis

The table below contains some constants you will need to complete the analysis described below

Variables	Description (Units)
Q	Flow Rate (ft^3/s)
\dot{m}	Mass Flow Rate (lb/s)
ρ	Density of Water (lb/ft^3)
A	Nozzle Area ($0.000845 ft^2$)
u	Velocity at Nozzle Exit (ft/s)
u_0	Velocity Deflected from Vane (ft/s)
s	Vertical Distance ($0.115 ft$)
F	Force on Vane ($lb.ft/s^2$)
M	Mass of Jockey Weight ($1.32 lb$)
g	Gravity (ft/s^2)
y	Horizontal Distance (ft)
a	Point of Application of F from Pivot ($0.50 ft$)

Flat Plate

1. Convert flow rate to mass flow rate, \dot{m}
2. Calculate the flow velocity, u
3. Calculate u_0
4. Calculate rate of delivery of momentum, $\dot{m}u_0$
5. Calculate force, F
6. Create a plot of force, F (y-axis) vs. rate of delivery of momentum, $\dot{m}u_0$ (x-axis)
 - Create a trendline for the data points
 - The slope of the trendline indicates the coefficient for theoretical force

Trial	Discharge (Q)	Mass Flow (\dot{m})	Distance (y)	Velocity (u)	Velocity (u_0)	Momentum Flux ($\dot{m}u_0$)	Force (F)
1							
2							
3							
4							
5							

Hemispherical Cup

1. Convert flow rate to mass flow rate, \dot{m}
2. Calculate the flow velocity, u
3. Calculate u_0
4. Calculate rate of delivery of momentum, $\dot{m}u_0$
5. Calculate force, F
6. Create a plot of force, F (y-axis) vs. rate of delivery of momentum, $\dot{m}u_0$ (x-axis)
 - Create a trendline for the data points
 - The slope of the trendline indicates the coefficient for theoretical force

Trial	Discharge (Q)	Mass Flow (\dot{m})	Distance (y)	Velocity (u)	Velocity (u_0)	Momentum Flux ($\dot{m}u_0$)	Force (F)
1							
2							
3							
4							
5							

Discussion/Interpretation

1. What percentage of velocity compared to that exiting the nozzle is lost as the jet deflects from the Flat Plate Vane?
2. What percentage of velocity compared to that exiting the nozzle is lost as the jet deflects from the Hemispherical Cup Vane?
3. Compare the coefficient of theoretical force found using the plot in Calculations 6b with the given, typical values for the coefficient of theoretical force on a Flat Plate Vane

- Compare your values for each trial using difference or some other mathematical comparison.
4. Compare the coefficient of theoretical force found using the plot in Calculations 6b with the given, typical values for the coefficient of theoretical force on a Hemispherical Cup Vane
- Compare your values for each trial using percentage difference or some other mathematical comparison
5. If an Angled Plate Vane, with an angle of 30° , was examined using Impact Jets, would you expect the coefficient of theoretical force to be higher or lower than the coefficients of theoretical force for Flat Plate and Hemispherical Cup Vanes?
- Justify your answer using data found in this experiment
6. Which source of error do you think is most significant in your experiment, and why?
7. Give an example of a real-life application of the principles and theories explored in this lab

Data Records

Record your readings in the tables already provided.