

# Laboratory 2 Sample Protocol

## Contents

- **Objective**
- **Materials and Equipment**
- **Procedure**
- **Data Analysis**
- **Deliverables**

## Objective

To investigate fluid statics by measuring buoyancy forces and hydrostatic thrust, and to validate theoretical principles using experimental data.

---

## Materials and Equipment

1. Quadrant balance apparatus
  2. Graduated cylinder (500 mL or larger)
  3. Thermometer
  4. Objects for buoyancy testing (rocks, composites, wood samples)
  5. Weighing scale ( $\pm 0.01$  g precision)
  6. Water (at ambient temperature)
  7. Ruler or measuring tape
  8. Transfer pipette
  9. Weight hangers and standard masses
-

# Procedure

## Part 1: Displacement Volumes and Buoyancy

### 1. Prepare the Setup:

- Measure and record the temperature of the water.
- Fill the graduated cylinder with a sufficient amount of water. Record the initial volume level, ( $V_{\text{initial}}$ ).

### 2. Measure Object Data:

- Weigh the first object (e.g., Rock-1) and record its mass.
- Gently submerge the object in the water and record the new volume level, ( $V_{\text{final}}$ ).
- Calculate the displaced volume, ( $\Delta V = V_{\text{final}} - V_{\text{initial}}$ ).
- Repeat the procedure three times for each object (Rock-2, Composite-1, etc.).

### 3. Repeat Measurements:

- Repeat the displacement experiment for all six objects.
- Ensure all measurements are consistent and record the data in a table.

## Part 2: Hydrostatic Forces and Center of Pressure

### 1. Prepare the Apparatus:

- Measure and record the water temperature.
- Verify both tanks in the quadrant balance are empty. Trim the assembly to ensure the submerged plane is vertical.

### 2. Partial Submersion:

- Add water into the trim tank to bring the balance to the 0 position. Add weights as needed to stabilize the apparatus.
- Gradually add water to the quadrant tank until the apparatus is level again. Record the water depth ( $h$ ) and the free surface width ( $b$ ).
- Repeat the procedure for at least three trials with varying weights.

### 3. Full Submersion:

- Fully submerge the plane surface by incrementally increasing weights and adding water to balance the apparatus.
- Record  $(h)$ ,  $(b)$ , and the applied masses for at least three trials.

## Data Analysis

### 1. Displacement Volumes and Buoyancy:

- Calculate the buoyancy force for each object using:  $[F_B = \rho_{\text{water}} \cdot \Delta V \cdot g]$
- Compare calculated object volumes with measurements from the displacement method.

### 2. Hydrostatic Forces:

- Calculate moments,  $(M)$ , using the formula:  $[M = W \cdot \left(\frac{3b}{8}\right) \cdot h]$
- Plot  $(M)$  vs.  $(h)$  for fully submerged data. Fit a straight line and compute  $(R^2)$ .
- Use the slope of the line to calculate the specific weight of water and compare it to literature values.

### 3. Partially Submerged Data:

- Plot:  $[M + \frac{\gamma_w W R_2^2 h}{2} \quad \text{vs.} \quad h^3]$
- Evaluate the fit using  $(R^2)$ .

## Deliverables

1. Completed data tables for Part 1 and Part 2.
2. Plots and calculations demonstrating experimental results.
3. A step-by-step experimental protocol with annotations for improvements.
4. Discussion addressing:
  - Archimedes' principle and its application.
  - Comparison of measured and theoretical buoyancy forces.
  - Analysis of hydrostatic forces and center of pressure.

